

Atlanta BioMedical Model 4100 Syringe Infusion Pump, Typical Spare Parts:




P/N 001-0022: Transducer Housing Seal



P/N 001-0040: Retainer Cluster Gear



P/N 001-0044: Seal/Cup



P/N 001-0050: Seal/Case (Rectangular)




P/N 001-0113: Seal/Cover AC Inlet




(With Protective Covers)

(The lens itself)


P/N 001-0116: Lens



P/N 001-00117: Keypad



P/N 001-0118: Overlay (White)



P/N 001-0122: Enteral Label Right (Purple) and 001-0123: Enteral Label Left (Purple)




P/N 001-0124: Enteral Label Right (Orange) and 001-0125: Enteral Label Left (Orange)



P/N 001-0142: Screw Battery Case (also see 111-0701 Seal)




P/N 001-0190: Overlay Orange



P/N 001-0191: Overlay (Purple)




P/N 001-0192: Medela Orange Overlay



P/N 001-0193: Spacer (LCD Lens)



P/N 001-0413: Transducer Seal




P/N 001-0431: Rubber Foot



P/N 001-0432: Ground Ring



P/N 001-0450: Clutch Lever



P/N 001-0490: Seal/Case (Circular)



P/N 001-5104: Battery Pack



P/N 001-5110: Main PCB Assembly



P/N 001-5111: Top PCB Assembly




P/N 001-5112: Track/Transducer Assembly



P/N 001-5114: Force Sensor Assembly




P/N 001-5116: Motor Assembly (use with S/N 5920 and below)



P/N 001-5117: LCD



P/N 001-5202: Potentiometer Assembly



P/N 001-5206: Top Cable Assembly



P/N 001-5209: AC Inlet Assembly



P/N 001-5211: Speaker Assembly




P/N 001-5217: Programming Cable (Cloning)



P/N 001-5218: Plunger-Retainer Assembly



P/N 001-5219: Track with Clutch



P/N 001-5223: Motor Assembly (use with 001-5224)

(use with S/N 5921 and above)



P/N 001-5224: Bottom with Saddle





P/N 001-5990: Pole Clamp Rotational



P/N 001-5991: Pump Handle



P/N 001-6010: Syringe Flange Retainer



P/N 001-6062: Cover, Transducer Housing



P/N 001-9200 and 001-9201, Inserts Set (Shipping)




(Side View)



(Bottom View)


P/N 001-9202: Box (Shipping)



P/N 111-0700: Seal/ Battery Case and 111-701: Seal/Battery Thumb Screw




P/N 111-0750: Seal/Plug Comm. Port




P/N 200-0061: Line Cord 6 ft (US Hosp Grade)



P/N 200-0062: Line Cord (220V)



P/N 200-0063: Line Cord 8 ft (US Hosp Grade)



P/N 200-0064: Line Cord Locking 10 ft (US Hosp Grade)



P/N 203-0104: Fuse .125A (120V unit)



P/N 203-0105: Fuse .080A (240V unit)



P/N 221-0401: Opto Sensor/Motor